




MARCH3 Antibody

Product Code	CSB-PA006853
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UD3
Immunogen	Synthesized peptide derived from the Internal region of Human 40605.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MARCH3; RNF173; E3 ubiquitin-protein ligase MARCH3; Membrane-associated RING finger protein 3; Membrane-associated RING-CH protein III; MARCH-III; RING finger protein 173
Immunogen Species	Homo sapiens (Human)
Target Names	MARCH3

Image



Western Blot analysis of K562 cells using MARCH3 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.