



RPL39L Antibody

Product Code	CSB-PA006881
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EH5
Immunogen	Synthesized peptide derived from the Internal region of Human Ribosomal Protein L39L.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RPL39L; RPL39L1; 60S ribosomal protein L39-like; 60S ribosomal protein L39-2
Immunogen Species	Homo sapiens (Human)
Target Names	RPL39L