




OR51F2 Antibody

Product Code	CSB-PA006909
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NH61
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 51F2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR51F2; Olfactory receptor 51F2; Olfactory receptor OR11-23
Immunogen Species	Homo sapiens (Human)
Target Names	OR51F2

Image



Western Blot analysis of K562 cells using Olfactory receptor 51F2 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.