




OR8I2 Antibody

Product Code	CSB-PA006932
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N0Y5
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 8I2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR8I2; Olfactory receptor 8I2; Olfactory receptor OR11-170
Immunogen Species	Homo sapiens (Human)
Target Names	OR8I2

Image



Western Blot analysis of HeLa cells using Olfactory receptor 8I2 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.