



CSNK1A1L/CSNK1A1 Antibody

Product Code	CSB-PA006977
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N752/P48729
Immunogen	Synthesized peptide derived from the Internal region of Human Casein Kinase I α 1/1L.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CSNK1A1L; Casein kinase I isoform alpha-like; CKI-alpha-like; CK1; CSNK1A1; Casein kinase I isoform alpha; CKI-alpha; CK1
Immunogen Species	Homo sapiens (Human)
Target Names	CSNK1A1L/CSNK1A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.