



CCR4 Antibody

Product Code	CSB-PA006992
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51679
Immunogen	Synthesized peptide derived from the Internal region of Human CKR-4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CCR4; CMKBR4; C-C chemokine receptor type 4; C-C CKR-4; CC-CKR-4; CCR-4; CCR4; K5-5; CD antigen CD194
Immunogen Species	Homo sapiens (Human)
Target Names	CCR4