




CCR5 Antibody

Product Code	CSB-PA006994
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51681
Immunogen	Synthesized peptide derived from Human CKR-5 around the non-phosphorylation site of S349.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CCR5; CMKBR5; C-C chemokine receptor type 5; C-C CKR-5; CC-CKR-5; CCR-5; CCR5; CHEMR13; HIV-1 fusion coreceptor; CD antigen CD195
Immunogen Species	Homo sapiens (Human)
Target Names	CCR5

Image



Western Blot analysis of COS7 cells using CCR-5 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.