




USP43 Antibody

Product Code	CSB-PA007025
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q70EL4
Immunogen	Synthesized peptide derived from the N-terminal region of Human USP43.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	USP43; Ubiquitin carboxyl-terminal hydrolase 43; Deubiquitinating enzyme 43; Ubiquitin thioesterase 43; Ubiquitin-specific-processing protease 43
Immunogen Species	Homo sapiens (Human)
Target Names	USP43

Image



Western Blot analysis of JK cells using USP43 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.