



Phospho-HSPB6 (S16) Antibody

Product Code	CSB-PA007033
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14558
Immunogen	Synthesized peptide derived from Human HSP20 around the phosphorylation site of S16.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HSPB6; Heat shock protein beta-6; HspB6; Heat shock 20 kDa-like protein p20
Immunogen Species	Homo sapiens (Human)
Target Names	HSPB6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.