



HSPB6 Antibody

Product Code	CSB-PA007034
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14558
Immunogen	Synthesized peptide derived from Human HSP20 around the non-phosphorylation site of S16.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HSPB6; Heat shock protein beta-6; HspB6; Heat shock 20 kDa-like protein p20
Immunogen Species	Homo sapiens (Human)
Target Names	HSPB6