




OR2AJ1 Antibody

Product Code	CSB-PA007058
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGZ0
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 2AJ1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR2AJ1; OR2AJ1P; Olfactory receptor 2AJ1
Immunogen Species	Homo sapiens (Human)
Target Names	OR2AJ1

Image



Western Blot analysis of Jurkat cells using Olfactory receptor 2AJ1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.