





COL4A2 Antibody

Product Code	CSB-PA007068
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08572
Immunogen	Synthesized peptide derived from the N-terminal region of Human COL4A2.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	COL4A2; Collagen alpha-2(IV) chain
Immunogen Species	Homo sapiens (Human)
Target Names	COL4A2

Image



Western Blot analysis of SH-SY5Y cells using COL4A2 Polyclonal Antibody



Western Blot analysis of SH-SY5Y 293T
NIH-3T3 COLO cells using COL4A2 Polyclonal
Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.