



GAB4/GAB2 Antibody

Product Code	CSB-PA007075
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2WGN9/Q9UQC2
Immunogen	Synthesized peptide derived from the N-terminal region of Human Gab 2/4.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GAB4; GRB2-associated-binding protein 4; GRB2-associated binder 2-like; GAB2-like; GRB2-associated binder 4; GRB2-associated-binding protein 2-like; Growth factor receptor bound protein 2-associated protein 4; GAB2; KIAA0571; GRB2-associate
Immunogen Species	Homo sapiens (Human)
Target Names	GAB4/GAB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.