



COL9A3 Antibody

Product Code	CSB-PA007102
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14050
Immunogen	Synthesized peptide derived from the Internal region of Human COL9A3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	COL9A3; Collagen alpha-3(IX) chain
Immunogen Species	Homo sapiens (Human)
Target Names	COL9A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.