



GLYCTK Antibody

Product Code	CSB-PA007152
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IVS8
Immunogen	Synthesized peptide derived from the Internal region of Human Glycerate Kinase.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GLYCTK; HBEBP4; LP5910; Glycerate kinase; HBeAg-binding protein 4
Immunogen Species	Homo sapiens (Human)
Target Names	GLYCTK