



H1FOO Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA007154 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8IZA3 |
| Immunogen | Synthesized peptide derived from Human Histone H1 around the non-acetylation site of K25. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | H1FOO; H1OO; OSH1; Histone H1oo; Oocyte-specific histone H1; Oocyte-specific linker histone H1; osH1 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | H1FOO |