





COX10 Antibody

Product Code	CSB-PA007188
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12887
Immunogen	Synthesized peptide derived from the Internal region of Human COX10.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	COX10; Protoheme IX farnesyltransferase; mitochondrial; Heme O synthase
Immunogen Species	Homo sapiens (Human)
Target Names	COX10

Image



Western blot analysis of KB lysis using COX10 antibody.



Western blot analysis of mouse-kidney mouse-brain HeLa KB 293T lysis using COX10 antibody.



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.