



Phospho-CLDN7 (Y210) Antibody

Product Code	CSB-PA007233
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95471
Immunogen	Synthesized peptide derived from Human Claudin-7 around the phosphorylation site of Y210.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CLDN7; CEPTRL2; CPETRL2; Claudin-7; CLDN-7
Immunogen Species	Homo sapiens (Human)
Target Names	CLDN7