




APOOL Antibody

Product Code	CSB-PA007280
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXV4
Immunogen	Synthesized peptide derived from the C-terminal region of Human ApoOL.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	APOOL; CXorf33; FAM121A; Apolipoprotein O-like; Protein FAM121A
Immunogen Species	Homo sapiens (Human)
Target Names	APOOL

Image



Western blot analysis of SH-SY5Y 293T 3T3 lysis using ApoOL antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.