



PARP4 Antibody

Product Code	CSB-PA007321
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJK3
Immunogen	Synthesized peptide derived from the Internal region of Human PARP-4.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PARP4; ADPRTL1; KIAA0177; PARPL; Poly [ADP-ribose] polymerase 4; PARP-4; 193 kDa vault protein; ADP-ribosyltransferase diphtheria toxin-like 4; ARTD4; PARP-related/lalpal-related H5/proline-rich; PH5P; Vault poly(ADP-ribose) polymerase; VP
Immunogen Species	Homo sapiens (Human)
Target Names	PARP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.