





# CSF1R Antibody

<b>Product Code</b>	CSB-PA007334
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P07333
<b>Immunogen</b>	Synthesized peptide derived from Human c-Fms around the non-phosphorylation site of Y561.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB, IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	CSF1R; FMS; Macrophage colony-stimulating factor 1 receptor; CSF-1 receptor; CSF-1-R; CSF-1R; M-CSF-R; Proto-oncogene c-Fms; CD antigen CD115
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	CSF1R

## Image



Western Blot analysis of COLO cells using c-Fms Polyclonal Antibody



Western Blot analysis of HEPG2 using c-Fms Polyclonal Antibody

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.