




Phospho-CSK (S364) Antibody

Product Code	CSB-PA007348
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41240
Immunogen	Synthesized peptide derived from Human Csk around the phosphorylation site of S364.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CSK; Tyrosine-protein kinase CSK; C-Src kinase; Protein-tyrosine kinase CYL
Immunogen Species	Homo sapiens (Human)
Target Names	CSK

Image



Western Blot analysis of VEC cells using Phospho-Csk (S364) Polyclonal Antibody