




CSK Antibody

Product Code	CSB-PA007350
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41240
Immunogen	Synthesized peptide derived from the N-terminal region of Human Csk.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CSK; Tyrosine-protein kinase CSK; C-Src kinase; Protein-tyrosine kinase CYL
Immunogen Species	Homo sapiens (Human)
Target Names	CSK

Image



Western Blot analysis of Jurkat cells using Csk Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.