




Phospho-CSNK1A1/CSNK1A1L (Y294) Antibody

Product Code	CSB-PA007356
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48729/Q8N752
Immunogen	Synthesized peptide derived from Human Casein Kinase I α around the phosphorylation site of Y294.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CSNK1A1; Casein kinase I isoform alpha; CKI-alpha; CK1; CSNK1A1L; Casein kinase I isoform alpha-like; CKI-alpha-like; CK1
Immunogen Species	Homo sapiens (Human)
Target Names	CSNK1A1/CSNK1A1L

Image



Western Blot analysis of HBE cells using Phospho-Casein Kinase I α (Y294) Polyclonal Antibody