




CSNK1D Antibody

Product Code	CSB-PA007359
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48730
Immunogen	Synthesized peptide derived from the Internal region of Human Casein Kinase Iδ.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CSNK1D; HCKID; Casein kinase I isoform delta; CKI-delta; CKId; Tau-protein kinase CSNK1D
Immunogen Species	Homo sapiens (Human)
Target Names	CSNK1D

Image



Western Blot analysis of brain cells using Casein Kinase Iδ Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.