




CSNK1E Antibody

Product Code	CSB-PA007361
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49674
Immunogen	Synthesized peptide derived from the Internal region of Human Casein Kinase Iε.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CSNK1E; Casein kinase I isoform epsilon; CKI-epsilon; CKIε
Immunogen Species	Homo sapiens (Human)
Target Names	CSNK1E

Image



Western Blot analysis of LOVO cells using Casein Kinase Iε Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.