




VCAN Antibody

Product Code	CSB-PA007374
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13611
Immunogen	Synthesized peptide derived from the N-terminal region of Human Versican.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	VCAN; CSPG2; Versican core protein; Chondroitin sulfate proteoglycan core protein 2; Chondroitin sulfate proteoglycan 2; Glial hyaluronate-binding protein; GHAP; Large fibroblast proteoglycan; PG-M
Immunogen Species	Homo sapiens (Human)
Target Names	VCAN

Image



Western Blot analysis of mouse-brain cells using Versican Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.