




ADRA1B Antibody

Product Code	CSB-PA007387
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35368
Immunogen	Synthesized peptide derived from the C-terminal region of Human AR- α 1B.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ADRA1B; Alpha-1B adrenergic receptor; Alpha-1B adrenoreceptor; Alpha-1B adrenoceptor
Immunogen Species	Homo sapiens (Human)
Target Names	ADRA1B

Image



Western Blot analysis of JK cells using AR- α 1B Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.