




ADRA1A Antibody

Product Code	CSB-PA007406
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35348
Immunogen	Synthesized peptide derived from the C-terminal region of Human AR- α 1A.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ADRA1A; ADRA1C; Alpha-1A adrenergic receptor; Alpha-1A adrenoreceptor; Alpha-1A adrenoceptor; Alpha-1C adrenergic receptor; Alpha-adrenergic receptor 1c
Immunogen Species	Homo sapiens (Human)
Target Names	ADRA1A

Image



Western Blot analysis of IOVO cells using AR- α 1A Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.