




CTBP1 Antibody

Product Code	CSB-PA007413
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13363
Immunogen	Synthesized peptide derived from Human CtBP1 around the non-phosphorylation site of S422.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CTBP1; CTBP; C-terminal-binding protein 1; CtBP1
Immunogen Species	Homo sapiens (Human)
Target Names	CTBP1

Image



Western Blot analysis of Jurkat cells using CtBP1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.