



Phospho-CTNNA1 (S641) Antibody

Product Code	CSB-PA007417
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35221
Immunogen	Synthesized peptide derived from Human Catenin- α 1 around the phosphorylation site of S641.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CTNNA1; Catenin alpha-1; Alpha E-catenin; Cadherin-associated protein; Renal carcinoma antigen NY-REN-13
Immunogen Species	Homo sapiens (Human)
Target Names	CTNNA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.