



# Phospho-CTNNB1 (Y489) Antibody

Product Code	CSB-PA007422
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35222
Immunogen	Synthesized peptide derived from Human Catenin- $\beta$ around the phosphorylation site of Y489.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CTNNB1; CTNNB; OK/SW-cl.35; Catenin beta-1; Beta-catenin
Immunogen Species	Homo sapiens (Human)
Target Names	CTNNB1