




ADRA2A Antibody

Product Code	CSB-PA007425
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08913
Immunogen	Synthesized peptide derived from the Internal region of Human AR- α 2A.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ADRA2A; ADRA2R; ADRAR; Alpha-2A adrenergic receptor; Alpha-2 adrenergic receptor subtype C10; Alpha-2A adrenoreceptor; Alpha-2A adrenoceptor; Alpha-2AAR
Immunogen Species	Homo sapiens (Human)
Target Names	ADRA2A

Image



Western Blot analysis of JK cells using AR- α 2A Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.