





FAM117B Antibody

Product Code	CSB-PA007427
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P1L5
Immunogen	Synthesized peptide derived from the C-terminal region of Human ALS2CR13.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FAM117B; ALS2CR13; Protein FAM117B; Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 13 protein
Immunogen Species	Homo sapiens (Human)
Target Names	FAM117B

Image



Western Blot analysis of HELA cells using ALS2CR13 Polyclonal Antibody



Western Blot analysis of HELA cells using ALS2CR13 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.