



PPP1R1C Antibody

Product Code	CSB-PA007453
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVI7
Immunogen	Synthesized peptide derived from the N-terminal region of Human PPP1R1C.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PPP1R1C; Protein phosphatase 1 regulatory subunit 1C; Inhibitor-5 of protein phosphatase 1; IPP5
Immunogen Species	Homo sapiens (Human)
Target Names	PPP1R1C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.