



GPBAR1 Antibody

Product Code	CSB-PA007456
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDU6
Immunogen	Synthesized peptide derived from the N-terminal region of Human TGR5.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; WB:1:500-2000, ELISA:2000-20000, IHC:1:50-300
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GPBAR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.