




PPP4R2 Antibody

Product Code	CSB-PA007474
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NY27
Immunogen	Synthesized peptide derived from the Internal region of Human PP4R2.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PPP4R2; SBBI57; Serine/threonine-protein phosphatase 4 regulatory subunit 2
Immunogen Species	Homo sapiens (Human)
Target Names	PPP4R2

Image



Western Blot analysis of K562 cells using PPP4R2 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.