



CNTN4 Antibody

Product Code	CSB-PA007480
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWW2
Immunogen	Synthesized peptide derived from the Internal region of Human Contactin 4.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CNTN4; Contactin-4; Brain-derived immunoglobulin superfamily protein 2; BIG-2
Immunogen Species	Homo sapiens (Human)
Target Names	CNTN4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.