



CNTN4 Antibody

Product Code	CSB-PA007480
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8I WV2
Immunogen	Synthesized peptide derived from the Internal region of Human Contactin 4.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CNTN4; Contactin-4; Brain-derived immunoglobulin superfamily protein 2; BIG-2
Immunogen Species	Homo sapiens (Human)
Target Names	CNTN4