




CYP2A6V2 Antibody

Product Code	CSB-PA007527
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13120
Immunogen	Synthesized peptide derived from the Internal region of Human CYP2A6V2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	Cytochrome P450
Immunogen Species	Homo sapiens (Human)
Target Names	CYP2A6V2

Image



Western Blot analysis of L929 cells using CYP2A6V2 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.