





ADRB3 Antibody

Product Code	CSB-PA007545
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13945
Immunogen	Synthesized peptide derived from the Internal region of Human AR-β3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ADRB3; ADRB3R; B3AR; Beta-3 adrenergic receptor; Beta-3 adrenoreceptor; Beta-3 adrenoceptor
Immunogen Species	Homo sapiens (Human)
Target Names	ADRB3

Image



Western Blot analysis of K562 cells using AR-β3 Polyclonal Antibody



Western Blot analysis of 22RV1 cells using AR-β3 Polyclonal Antibody