



ADRBK1 Antibody

Product Code	CSB-PA007559
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25098
Immunogen	Synthesized peptide derived from Human GRK 2 around the non-phosphorylation site of S29.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ADRBK1; BARK; BARK1; GRK2; Beta-adrenergic receptor kinase 1; Beta-ARK-1; G-protein coupled receptor kinase 2
Immunogen Species	Homo sapiens (Human)
Target Names	GRK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.