




CYP2E1 Antibody

Product Code	CSB-PA007577
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05181
Immunogen	Synthesized peptide derived from the C-terminal region of Human CYP2E1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CYP2E1; CYP2E; Cytochrome P450 2E1; 4-nitrophenol 2-hydroxylase; CYP1IE1; Cytochrome P450-J
Immunogen Species	Homo sapiens (Human)
Target Names	CYP2E1

Image



Western Blot analysis of KB cells using CYP2E1 Polyclonal Antibody