




NKX6-3 Antibody

Product Code	CSB-PA007583
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NJ46
Immunogen	Synthesized peptide derived from the Internal region of Human Nkx-6.3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Monkey
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NKX6-3; Homeobox protein Nkx-6.3
Immunogen Species	Homo sapiens (Human)
Target Names	NKX6-3

Image



Western Blot analysis of K562 cells using NKx-6.3 Polyclonal Antibody