




ZDHHC15 Antibody

Product Code	CSB-PA007602
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96MV8
Immunogen	Synthesized peptide derived from the C-terminal region of Human DHHC-15.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZDHHC15; Palmitoyltransferase ZDHHC15; Zinc finger DHHC domain-containing protein 15; DHHC-15
Immunogen Species	Homo sapiens (Human)
Target Names	ZDHHC15

Image



Western Blot analysis of Jurkat cells using DHHC-15 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.