




DGKB Antibody

Product Code	CSB-PA007623
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6T7
Immunogen	Synthesized peptide derived from the C-terminal region of Human DGK-β.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DGKB; DAGK2; KIAA0718; Diacylglycerol kinase beta; DAG kinase beta; 90 kDa diacylglycerol kinase; Diglyceride kinase beta; DGK-beta
Immunogen Species	Homo sapiens (Human)
Target Names	DGKB

Image



Western Blot analysis of HuvEC cells using DGK-β Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.