




DGKQ Antibody

Product Code	CSB-PA007638
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52824
Immunogen	Synthesized peptide derived from the C-terminal region of Human DGK-θ.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DGKQ; DAGK4; Diacylglycerol kinase theta; DAG kinase theta; Diglyceride kinase theta; DGK-theta
Immunogen Species	Homo sapiens (Human)
Target Names	DGKQ

Image



Western Blot analysis of HT-29 cells using DGK-θ Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.