




DAPK1 Antibody

Product Code	CSB-PA007646
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53355
Immunogen	Synthesized peptide derived from Human DAPK1 around the non-phosphorylation site of S308.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DAPK1; DAPK; Death-associated protein kinase 1; DAP kinase 1
Immunogen Species	Homo sapiens (Human)
Target Names	DAPK1

Image



Western Blot analysis of MOUSE-BRAIN cells using DAPK1 Polyclonal Antibody