




Phospho-DAPK3 (T265) Antibody

Product Code	CSB-PA007652
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43293
Immunogen	Synthesized peptide derived from Human DAPK3 around the phosphorylation site of T265.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DAPK3; ZIPK; Death-associated protein kinase 3; DAP kinase 3; DAP-like kinase; Dlk; MYPT1 kinase; Zipper-interacting protein kinase; ZIP-kinase
Immunogen Species	Homo sapiens (Human)
Target Names	DAPK3

Image



Western Blot analysis of HuvEc cells using Phospho-DAPK3 (T265) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.