




TCEB3C Antibody

Product Code	CSB-PA007655
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NG57
Immunogen	Synthesized peptide derived from the C-terminal region of Human Elongin A3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TCEB3C; TCEB3L2; RNA polymerase II transcription factor SIII subunit A3; Elongin-A3; EloA3; Transcription elongation factor B polypeptide 3C
Immunogen Species	Homo sapiens (Human)
Target Names	ELOA3

Image



Western Blot analysis of COLO205 cells using Elongin A3 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.